

Blood Banking & Transfusion Medicine 101

ABO Typing, Antibody Screening, Identification & Crossmatching Why is This Important

Presented by:
Adam Norfolk, MBA, MLS(ASCP)^{CM}SBB^{CM}
Lead Technologist Blood Bank
Beth Israel Deaconess Medical Center
Boston, MA

www.aabb.org

1

Learning Objectives

After participating in this program you should be able to....

- Reveal the major ABO types and their prevalence in the donor and patient population
- Explain what front (forward) and back (reverse) types are and how to perform each
- Describe the purpose of the antibody screen
 - Explain why group O cells are always used in an antibody screen
 - Describe two cell versus three cell screens
- Explain the purpose of antibody identification
- Describe what a crossmatch is
 - Explain the differences between immediate spin, full serological (AHG) and electronic crossmatch




www.aabb.org

2

2

What is Blood?


- Blood has 4 major components:
 1. Red Blood Cells (RBCs)
 2. Plasma
 3. White Blood Cells (WBCs)
 4. Platelets



3

What is Blood?


- Blood Bank Tests:
 1. Red Blood Cells (RBCs)
 - Antigens
 2. Plasma
 - Where antibodies are stored



4

Antigens vs Antibodies

- Antigens if foreign to the body, can stimulate the immune system. Examples:
 - Bacteria
 - Virus
 - RBC
- A stimulated immune system may make antibodies to interact with the antigen.



5

What is an ABO?

- An *ABO* or *Blood Type* references the ABO blood group system, one of 36 blood group systems!
 - Hundreds of antigens on RBCs.
- The ABO blood group system is the most important blood group system.
- The four blood types are: O, A, B, AB



1. Fung M, Eder A, Spitalnik S, Westhoff C. AABB Technical Manual. 19th ed. Bethesda: AABB. 2017; 262 p.

6

What is an ABO?

	Type O	Type A	Type B	Type AB
Antigen		A	B	AB
Antibody	Anti-A Anti-B	Anti-B	Anti-A	No antibodies
Transfuse RBC ABO Group	O	A or O	B or O	AB or A or B or O
Transfuse Plasma ABO Group	O or A or B or AB	A or AB	B or AB	AB



7

ABO Prevalence (%) in the US Population

ABO Group	European Ethnicity	African Ethnicity
O	45	49
A	40	27
B	11	20
AB	4	4


AABB Technical Manual, 19 ed. Table 10-1 page 267



8

Blood Bank Testing

- Hold
- Type and Hold
 - “Type” or “ABO/Rh” or “Blood Type”
- Type and Screen
 - ABO & Rh Type
 - Antibody Screen
 - Antibody Identification if needed
- Crossmatch



9

ABO & Rh Type

- Front Type
 - Determines which antigens are on the RBCs
 - Reagent ABO antibody bind to patient RBC antigens.



10

ABO Front Type

- 1 Place 1 drop of anti-A in a clean, labeled test tube.
- 2 Place 1 drop of anti-B in a separate, clean, labeled tube.
- 3 To each tube, add 1 drop of a 2% to 5% suspension (in saline, serum, or plasma) of the red cells to be tested. Alternatively, the equivalent amount of red cells can be transferred to each tube with clean applicator sticks.
- 4 Gently mix the contents of the tubes; then centrifuge for the calibrated spin time.
- 5 Gently resuspend the cell buttons, and examine them for agglutination.
- 6 Read, interpret, and record the test results. Compare the red cell test results with those obtained in the serum or plasma tests.




www.aabb.org

11

11


Testing Principles-Agglutination




12

Macroscopic Agglutination

Negative



Positive




www.aabb.org

13


13

Microscopic Agglutination

Negative



Positive



www.aabb.org

14

14

ABO Backtype

- Back Type
 - Determines which antibodies are present in the patient
 - Patient's antibody binds to reagent red cell antigens.



15

ABO Back Type

- | | |
|---|--|
| 1 | Add 2 or 3 drops each of serum or plasma to two clean, labeled test tubes. |
| 2 | Add 1 drop of A ₁ reagent red cells to the tube labeled A ₁ . |
| 3 | Add 1 drop of B reagent red cells to the tube labeled B. |
| 4 | Gently mix the contents of the tubes; then centrifuge for the calibrated spin time. |
| 5 | Examine the serum overlying the red cell buttons for evidence of hemolysis. Gently resuspend the cell buttons, and examine them for agglutination. |
| 6 | Read, interpret, and record test results. Compare serum test results with those obtained in testing red cells. |




www.aabb.org

16

16

Step 1: Label Tubes




www.aabb.org

17

17

Step 1: Label Tubes



www.aabb.org

18

18

Step 2: Add Saline to Cell Suspension Tube




www.aabb.org

19

19

Step 3: Obtain Patient Plasma




www.aabb.org

20

20

Step 4: Add 2 Drops of Plasma to Backtype Tubes




www.aabb.org

21

21

Step 5: Prepare Patient Red Cell Suspension



www.aabb.org

22

22

Step 6: Add anti-A




www.aabb.org

23

23

Step 7: Add anti-B




www.aabb.org

24

24

Step 8: Add anti-D




www.aabb.org

25

25

Step 9: Add Reagent Red Cells to Back Type



www.aabb.org

26

26

Step 10: Add Reagent A1-cells




www.aabb.org

27

27

Step 11: Add Reagent B-Cells




www.aabb.org

28

28

Step 12: Add Patient Cells to Front Type




www.aabb.org

29

29

Step 13: Centrifuge




www.aabb.org

30

30

Step 14: Shake, Rattle, and Roll-Front Type




www.aabb.org

31

31

Step 14: Shake, Rattle, and Roll-Front Type




www.aabb.org

32

32

Step 14: Shake, Rattle, and Roll-Front Type




www.aabb.org

33

33

Step 14: Shake, Rattle, and Roll-Back Type



www.aabb.org

34

34

Step 14: Shake, Rattle, and Roll-Back Type




www.aabb.org

35

35

Step 14: Shake, Rattle, and Roll-D Type



www.aabb.org

36

36

ABO Grouping

Routine ABO Grouping

Reaction of Red Cells with Antisera (Red Cell Grouping)		Reaction of Serum with Reagent Red Cells (Serum Grouping)			Interpretation
Anti-A	Anti-B	A ₁ Cells	B Cells	O Cells	
0	0	+	+	0	O
+	0	0	+	0	A
0	+	+	0	0	B
+	+	0	0	0	AB



Fung M, Eder A, Spitalnik S, Westhoff C. AABB Technical Manual. 19th ed. Bethesda: AABB. 2017; Table 12-2.

www.aabb.org

37

37

Detecting RBC Antibodies



www.aabb.org

38

38

The Antibody Screen

- To detect clinically significant antibodies to non-ABO antigens.
- Reagent red cells must be ABO group O.
- FDA requires the following antigens be present for licensed reagent screens:
 - D,C,E,c,e; K,k; Fy^a,Fy^b; Jk^a,Jk^b; Le^a,Le^b; P1; M,N,S,s



www.aabb.org

39

39

The Antibody Screen

- Licensed reagent antibody screens have 2 or 3 cells.
 - 2 cell screen
 - Both cells are D-Positive
 - Variable expression of other antigens.
 - 3 cell screen
 - 2 cells are D-Positive and 1 is D-Negative
 - Better differentiation of antibody if detected.
 - Usually offer stronger expression of antigens.



www.aabb.org

40

40

The Antibody Screen

2-Cell

	RH						MNS				LU		P	Lewis		Kell		Duffy		Kidd		Gel Test	
	D	C	E	c	e	f	M	N	S	s	Lu ^a	Lu ^b	P ₁	Le ^a	Le ^b	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	IAT	
1	+	+	0	0	+	0	+	+	+	+	0	+	+	+	0	+	+	+	0	0	+	1+	
2	+	0	+	+	0	0	+	0	0	+	0	+	+	0	+	0	+	+	+	+	0	0	

3-Cell

	RH						MNS				LU		P	Lewis		Kell		Duffy		Kidd		LISS (Test Tube)		
	D	C	E	c	e	f	M	N	S	s	Lu ^a	Lu ^b	P ₁	Le ^a	Le ^b	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	IS	37 C	IAT
1	+	+	0	0	+	0	+	+	+	+	0	+	+	+	0	+	+	+	0	0	+	0	0	2+
2	+	0	+	+	0	0	+	0	+	0	0	+	+	0	+	0	+	+	+	+	0	0	0	1+
3	0	0	0	+	+	+	0	+	0	+	0	+	0	0	+	0	+	0	+	+	+	0	0	0√



www.aabb.org

Antibody Identification: Art or Science? A Case Study Approach
 By Janis R. Hamilton, MS, MT(ASCP)SBB; Susan T. Johnson, MSTM, MT(ASCP)SBB;
 Sally V. Rudmann, PhD, MT(ASCP)SBB
 Copyright © 2013 by AABB. All rights reserved.

41

41

Two-Cell Antibody Screen


1. Add 2 drops of patient plasma into 2 clean tubes labeled "I" and "II".
2. Add 1 drop of screening cell I into tube "I".
3. Add 1 drop of screening cell II into tube "II".
4. Add 2 drops of enhancement reagent (PEG). Mix well.
5. Incubate at 37C for 15 minutes.
6. Wash the red cells four times with saline, and completely decant the final wash.
7. Add anti-IgG to the dry red cell button. Mix well.
8. Centrifuge and observe for agglutination. Grade and record the results.
9. Confirm the validity of negative results by adding IgG-coated red cells.




www.aabb.org

42


42



43



44



45

Antibody Identification

- Performed if the screen is positive.
- A series of reagent cells, typically 10, 16, or 20 cells, with variable phenotypes for the major blood group systems.
 - Allows for patterns of reactivity to rule in and rule out antibodies.
 - Enhancement media may be added to reduce 37C incubation time from 60 minutes to 15 minutes just like the antibody screen.
- It's a puzzle!



www.aabb.org

46

46

Antibody ID Antigram

	RH						MNS				LU		P	Lewis		Kell		Duffy		Kidd		SP	
	D	C	E	c	e	f	M	N	S	s	Lu ^a	Lu ^b	P ₁	Le ^a	Le ^b	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	IAT	
1	+	+	0	0	+	0	+	+	+	+	0	+	+	0	+	0	+	0	+	+	+	+	2+
2	+	+	0	0	+	0	+	+	+	0	0	+	0	0	0	+	+	+	+	+	+	+	2+
3	+	0	+	+	0	0	0	+	+	+	0	+	0	0	0	+	+	+	+	0	+	+	2+
4	+	0	0	+	+	+	+	+	0	+	0	+	+	+	0	0	+	0	0	0	0	+	2+
5	0	+	0	+	+	+	+	+	0	+	0	+	+	+	0	0	+	+	0	+	0	+	0
6	0	0	+	+	+	+	+	+	+	+	0	+	+	0	+	0	+	0	+	+	+	+	0
7	0	0	0	+	+	+	+	+	0	+	0	+	+	+	0	+	+	0	+	0	+	+	0
8	0	0	0	+	+	+	+	+	+	+	0	+	0	0	+	0	+	+	0	+	0	+	0
9	0	0	0	+	+	+	+	0	+	0	0	+	+	+	0	0	+	0	+	+	+	+	0
10	0	0	0	+	+	+	0	+	0	+	0	+	+	0	+	0	+	+	0	0	0	+	0
11	+	+	0	0	+	0	+	+	+	+	0	+	0	+	0	+	+	0	+	+	+	+	2+



www.aabb.org

Antibody Identification: Art or Science? A Case Study Approach
 By Janis R. Hamilton, MS, MT(ASCP)SBB; Susan T. Johnson, MSTM, MT(ASCP)SBB;
 Sally V. Rudmann, PhD, MT(ASCP)SBB
 Copyright © 2013 by AABB. All rights reserved.

47

47

Crossmatching Blood

How do we detect incompatibilities?



www.aabb.org

48

48

Crossmatch

- Immediate-Spin
 - Serological method to detect ABO incompatibility
- Computer/Electronic
 - Alternative way to verify ABO compatibility
- Antiglobulin
 - Serological method to detect RBC compatibility to non-ABO clinically significant antibodies.



www.aabb.org

49

49

Immediate-Spin Crossmatch


- Donor red cells and patient plasma are mixed, spun, and read.
- Can be used as the sole crossmatch method only if the recipient has no present or previously detected, non-ABO clinically significant antibodies.
- Failure to follow procedure can result in failure to detect ABO-incompatible RBCs.



www.aabb.org

50

50



51

Computer/Electronic Crossmatch

- Can be used if:
 - System is validated to ensure only ABO-compatible RBCs are given.
 - Two determinations of the patient's ABO group are made
 - The system contains the DIN, component name, ABO & Rh group, the confirmed donor unit's ABO group, two unique patient identifiers, the patient's ABO, Rh, ABS result and compatibility interpretation.
 - A method to verify correct entry of data before release of blood.
 - System contains logic to alert users of discrepancies or incompatibility.

The AABB logo and website are at the bottom left, and the slide number '52' is at the bottom right.

52

Antiglobulin Crossmatch

- Used for patient with historical or currently detected non-ABO clinically significant antibodies.
- RBC must lack relevant antigens the antibodies are directed to.
- Crossmatch mixing donor RBC and patient plasma, a 37C incubation, followed by the AHG test.
 - Enhancement media may be added to reduce 37C incubation time from 60 minutes to 15minutes just like the antibody screen and antibody ID!
- Not a substitute for the IS crossmatch.




www.aabb.org

53

53

Antiglobulin Crossmatch

- PEG
- Patient Plasma
- Red Cell Unit
- Anti-IgG
- Check Cells



www.aabb.org

54

54

Why Is This Important?

- Transfusing blood is the most common procedure performed in the USA.
- Support fetal development, surgeries, cancer treatments, CMO patient's.
- National Security
- Failure to provide safe transfusions can be fatal.



www.aabb.org

55

55

Questions?

Contact
AABB eLearning Team
eLearning@aabb.org



www.aabb.org

Program is copyright 2020 by AABB with All Rights Reserved

56

56